


High Resolution Automatic Viscometer Distinguishes Injectability of Candidate Proteins at an Early Stage of Development

There has been a strong market need for instruments that can screen candidate protein therapeutics at an early stage of development, thereby allowing identification of risks that might otherwise not be apparent until later stages of development. Since candidate compounds are usually available in small quantity, the ability to test small samples is critically important.

VROC initium™ is designed precisely to meet these stringent requirements, while providing exceptional resolution in viscosity measurement in order to distinguish viscosity differences as small as 0.020 mPa-s. Additionally, automation allows high throughput viscosity measurement of many candidate protein therapeutics on an unattended basis. The table and graph below demonstrate the excellent resolution and repeatability achieved; small differences in viscosity due to concentration differences in a dilute regime are accurately measured.

BSA			Bovine γ -Globulin		
Concentration, mg/mL	Viscosity, mPa-s	C_v	Concentration, mg/mL	Viscosity, mPa-s	C_v
10.34	0.960	0.0042	10.73	0.982	0.0053
19.93	1.005	0.0069	21.48	1.068	0.0066
28.61	1.050	0.0070	28.89	1.134	0.0133
35.64	1.088	0.0062	35.67	1.202	0.0053
47.64	1.158	0.0085	47.60	1.354	0.0074

C_v : Coefficient of variation defined as standard deviation divided by average value.



The test was conducted as follows:

Experiment:





1. 100 μL aliquots of protein solution samples are loaded into vial inserts, which are then placed into 11 mm vials.
2. The vials are placed into a 40 position rack on the VROC initium.
3. The Autosampler picks up 80 μL of sample and injects the sample into a test syringe for automated viscosity measurement.
4. The VROC® initium unit performs viscosity measurement at a shear rate of 9,581 1/s seven times for each sample, following preset protocols.

If you have questions or need more information about this product or other applications, please contact us:

Main Office — 1 925 866 3801
Information — info@RheoSense.com
Sales — Sales@RheoSense.com

