

microVISC

Chip	Flow Channel Depth (μm)	Viscosity Range (cP)
A02	20	0.2 - 4
A05	50	0.3 - 80
A10	100	2 - 600
A20	200	20 - 5,000
A30	300	70 - 17,000
B02	20	0.2 - 20
B05	50	1 - 300
B10	100	10 - 2,000
B20	200	70 - 16,000
B30	300	250 - 60,000
C05	50	6 - 1,400
C10	100	40 - 10,000
C20	200	350 - 80,000
C30	300	1,200 - 280,000
E20	200	3,000 - 650,000