



RheoSense

Simply Precise™

| Temperature Control Specifications (HVROC-T)



- ✓ Precise Temperature Control
- ✓ More Accurate Results
- ✓ Easy Keyboard Features

Portable *microVISCT*™ viscometer can be placed in the temperature control unit for precise controlled measurements.

Power	100-240 VAC; 50/60 Hz; 75W
Dimension	7.5 in x13in x15in
Temp Range	18-50 °C (64.4-122 °F)
Temp Stability	±.07 °C
Control Mechanism	Peltier Based Air Cooling/Heating
Controller	Fuzzy Logic PID
PC Interface	USB
Display	1.5x2.5 inch LCD
Temperature Sensor	100 ohm RTD



2420 Camino Ramon, Suite 240
San Ramon, CA 94583

P: (925) 866-3801 F: (925) 866-3804
RheoSense.com