

List of Consumables

Product Number	Description	Recommended Replacement	Price, \$ (USD)
INI-C-1004	100µl Syringe Assembly with an Injection Port	1000 uses	\$927.00
INI-C-1011	100µl Autosampler Syringe Assembly	1000 uses	\$1041.00
INI-C-1019	Reusable Injection Septum*	400 uses IP/700 uses WP	\$228.00
INI-C-1022	100 x 1.1 ml vials set with conical insert	As needed	\$103.00
INI-C-1023	100 x 1.1 mL vial caps	As needed	\$60.00
INI-C-1028	Aquet Pouch	As needed	\$12.00
H-300-4	MGVS20 Water based calibration fluid (2 cP) – 100 mL	As needed	\$150.00
H-300-5	MGVS60 Water based calibration fluid (6 cP) – 100 mL	As needed	\$150.00
H-300-6	MGVS100 Water based calibration fluid (10 cP) – 100 mL	As needed	\$150.00

*P/N INI-C-1019 has two different use counts that are based on the location that the septum is being used. If the septum is used in the injection port, the recommended use count is 400. If the septum is being used in the wash port, the recommended use count is 700.

List of Accessories & Replacement Parts

Product Number	Description	Price, \$ (USD)
INI-C-1008	100µl Manual Sample Injection Syringe	\$105.00
INI-C-1016	500 mL Media Bottle with Cap	\$208.00
INI-C-1017	250 mL Media Bottle with Cap	\$208.00
INI-C-1018	Two Ports Solvent Bottle Cap	\$151.00
INI-C-1020	40 Vial Rack	\$308.00
INI-C-1021	96 Well Plate Rack	\$371.00
INI-C-1024	peek NEEDLE (250µm ID) and fitting of the 100 µl test syringe (Qty 1)	\$39.00
INI-C-1025	PFA Outlet Tube for the Chip & Fitting (Qty 1)	\$39.00
INI-C-1026	6-32 1/16" Tube Fitting (Lab-Smith Fitting)	\$23.00
INI-C-1032	Wrench Tool Kit (1/8", 5/64", 3/32")	\$53.00
INI-C-1034	Kit Bottle Solvent Waste System – 4L	\$661.00
INI-C-1035	Solvent Waste HD Exhaust Filter	\$778.00
INI-P-1020	Power Adaptor	\$413.00
	Inlet Tube for Chip	
	Solvent Bottle Tube Set	

Service Contract Options

A service contract is the most economical method of keeping your Initium running at peak performance, while also providing peace of mind that all repairs of your instrument will be covered. At RheoSense, we offer service plans that cover the entire spectrum of customer needs. For more information about service plan options, please see the second page of this handout or contact a RheoSense Technical Specialist.

Product Number	Description	Price, \$ (USD)
INI-H-1003	Full Platinum Service Maintenance Plan	\$30,000
INI-H-1004	Full Gold Service Maintenance Plan	\$25,000
INI-H-1005	Full Silver Service Maintenance Plan	\$17,000

