


m-VROC II Chip Options for Regular Speed Configuration


WWW.RHEOSENSE.COM/PRODUCTS/VISCOMETERS/M-VROC-II



CHIP TYPE	SYRINGE SIZE	Viscosity Range, mPa-s		Shear Rate Range, 1/s	
		MIN	MAX	MIN	MAX
A05	100 μ L	0.04	80	40	22,000
B05	100 μ L	0.5	280	40	22,000
C05	100 μ L	3	1,500	40	22,000
C05	2.5 mL	1	70	900	558,000
C20	100 μ L	160	100,000	2	1,220
E02	100 μ L	1	80	1,300	143,000
E02	2.5 mL	1	6	17,000	2,020,000
E05	100 μ L	10	1,600	156	22,000
E05	2.5 mL	1	80	3,100	558,000

m-VROC II Chip Options for Fast Speed Configuration

WWW.RHEOSENSE.COM/PRODUCTS/VISCOMETERS/M-VROC-II



CHIP TYPE	SYRINGE SIZE	Viscosity Range, mPa-s		Shear Rate Range, 1/s	
		MIN	MAX	MIN	MAX
A05	100 μ L	0.04	80	80	44,000
B05	100 μ L	0.5	280	80	44,000
C05	100 μ L	3	1,500	80	44,000
C05	2.5 mL	1	70	1,800	558,000
C20	100 μ L	160	100,000	4	2,440
E02	100 μ L	1	80	1,300	287,000
E02	2.5 mL	0.4	6	17,000	5,000,000
E05	100 μ L	6	1,600	156	44,000
E05	2.5 mL	1	80	3,100	1,100,000