

Small Sample Viscosity Measurements

microVISC is the ultimate solution for your routine viscosity measurements




Small Sample Size

microVISC provides accurate and repeatable viscometry measurements using sample size as little as 100 microliters

Portable

Measuring at 14.5 x 21.6 x 7.1in and weighting under 2lbs, the microVISC can fit into a glove box for easy transport



Battery Operated

Battery operation allows you to run sample measurements anywhere! The rechargeable battery can run ~100 measurements between charges

Ease of Use

The microVISC is intuitively designed so that users can start operating within one minute



Simplified Data

Data collected can be automatically exported into a PDF report and/or Excel Spreadsheet formatting

Superior Accuracy

Accuracy up to 2% of reading for samples from low viscosity to high viscosity and repeatable measurement - 0.5% of reading

